



## Euroguard

Kirurgisk munnbind med visir

**Kirurgisk munnbind 160 med visir, Anti-gjenskinn og anti-dugg stripe.**

**Art. nr: 17203**

### Forpakning

25/boks, 100 per eske

Indre boks: 842942003738

Ytre boks: 10842942003735

**CE-merket**

Spesifikasjoner	
EN14683 klassifisering	Type IIR
BFE	> 99% (at 3.0 microns)
PFE	> 99% (at 0.1 microns)
Delta P	39.2 Pa
Væskemotstand	160 mm Hg (ASTM 1862)
Brennbarhet	Meets Class 1 (16 CFR 1610)
Dimensjoner	18.0 x 9.3 cm
Ytre lag	Polyspun bond polypropylene
Filter lag	Melt blown polypropylene
Innvendig fôr	Cellulose
Rør	Spunlace (Viscose/PET) top and bottom
Anti-dugg stripe:	Soft, open-celled foam
Anti-gjenskinn stripe:	Black PSB/PP
Visir	Anti-glare, anti-fog acrylic, with wrap-around wings
Konstruksjon	All seams ultrasonically welded – Total 6 layers
Nese-stykke	Aluminum alloy strip 14.0 x 0.46 cm
Øre-løkker	Spandex / polyester Rounded type 11.5 cm free length Ultrasonically welded to outer side

Testing: Performed at Nelson Labs, Salt Lake City (USA) in accordance with ASTM F2100-07 and EN14683 standards.